The cover features a white background with large, colorful curved shapes in blue, red, and green. A stylized logo above the text consists of a red dot, a green ring, a blue ring, and a red ring. The text 'LUXONUS' is in a bold, black, sans-serif font, and 'Technology guide' is in a smaller, black, sans-serif font below it. On the left, there is a black desk with a monitor displaying a technical diagram and a keyboard. In the center, there is a large, black, rectangular device with a small, colorful marker on top and a small, square, gold-colored feature on its surface.

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システム	Nagae, K. et al. Real-time 3D Photoacoustic Visualization System with a Wide Field of View for Imaging Human Limbs. F1000Research 7:1813 (2019).	
乳がん	Matsumoto, Y. et al. Visualising peripheral arterioles and venules through high-resolution and large-area photoacoustic imaging. Sci. Rep. 8 (2018).	
	Toi M, Asao Y, Matsumoto Y, et al.: Visualization of tumor-related blood vessels in human breast by photoacoustic imaging system with a hemispherical detector array. Sci Rep. 2017; 7: 41970.	
	Yamaga I et al., Vascular branching point counts using photoacoustic imaging in the superficial layer of the breast: A potential biomarker for breast cancer. Photoacoustics Volume 11, September 2018, Pages 6-13	
リンパ浮腫	Kajita H. and Kishi K., High-Resolution Imaging of Lymphatic Vessels with Photoacoustic Lymphangiography, Radiology 2019; 292:35	
	Kajita H. et al., Photoacoustic lymphangiography, J Surg Oncol. 2020; 121: 48-50	
	Suzuki Y. et al., Subcutaneous Lymphatic Vessels in the Lower Extremities: Comparison between Photoacoustic Lymphangiography and Near-Infrared Fluorescence Lymphangiography, Radiology VOL. 295, NO. 2	
遊離皮弁	Tsuge I. et al., Photoacoustic Tomography Shows the Branching Pattern of Anterolateral Thigh Perforators In Vivo, Plastic and Reconstructive Surgery: May 2018 - Volume 141 - Issue 5 - p 1288-1292	
	Tsuge I. et al., Preoperative vascular mapping for anterolateral thigh flap surgeries: A clinical trial of photoacoustic tomography imaging, Microsurgery, Volume40, Issue3, March 2020, Pages 324-330	
乾癬	Ishida Y., Photoacoustic imaging system visualizes restoration of peripheral oxygenation in psoriatic lesions, J Eur Acad Dermatol Venereol. 2018 Dec;32(12):e449-e451.	
手の血管	Matsumoto Y. et al., Label-free photoacoustic imaging of human palmar vessels: a structural morphological analysis, Scientific Reports volume 8, Article number: 786 (2018)	
	Saito S. et al., Digital artery deformation on movement of the proximal interphalangeal joint, Journal of Hand Surgery, Volume: 44 issue: 2, page(s): 187-195	

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